



GLIDESCOPE TITANIUM

REUSABLE VIDEO LARYNGOSCOPES



GlideScope® Titanium video laryngoscopes combine all-new blade designs and construction, making them our thinnest, toughest video laryngoscopes ever.

These new low-profile designs are lightweight and streamlined, offering improved maneuverability and working space for routine and difficult intubations.

The LoPro T3 and T4 combine the performance of a GlideScope blade with the strength of titanium. The MAC T3 and T4 provide the comfort of a Mac blade with the confidence of a GlideScope.

HIGH-PERFORMANCE FEATURES

- Choice of blade sizes and designs for a wide array of clinical demands
- Onboard 4-Step Technique video tutorial making the system easy to learn and use
- Real-time snapshot, video recording, and video playback features to help confirm tube placement and facilitate teaching
- Easy-to-clean blades – IPX8 compliant, with disinfection and sterilization options
- External monitor connection capability for education and observation

THE TITANIUM SYSTEM

- 4 reusable blades
- Low-profile designs for improved maneuverability and working space
- Titanium construction for lightweight durability
- New Mac-style blades for a wide range of clinical needs
- Digital, full-color camera and video monitor
- Reveal™ anti-fog technology, with a rapid heating profile to resist lens fogging
- Fully compatible with digital AVL blade configurations and monitors



GLIDESCOPE TITANIUM

REUSABLE VIDEO LARYNGOSCOPE

APPLICATIONS

Utilizing a laryngoscopic blade with an integrated camera connected to a video monitor, GlideScope® Titanium video laryngoscopes provide a clear view of the airway and vocal cords in procedures such as first-use intubations, normal or restricted oropharyngeal views/visualization, intubations with blood or other secretions in the airway, airway management in morbidly obese patients, preterm and neonatal intubations, as well as routine intubations.

SYSTEM COMPONENTS

- Reusable, titanium blade with digital camera
- Digital, full-color, rechargeable video monitor
- Video cable for connecting the blade to the monitor
- GlideRite® Rigid Stylet
- Choice of mobile stand or IV pole mount
- In-service USB flash drive including the Operations & Maintenance Manual
- 12-volt DC power adapter and power cable
- (Optional) HDMI-to-DVI cable for connecting to an external monitor
- Compatible AVL pediatric options available

VIDEO LARYNGOSCOPES & ADDITIONAL SYSTEM COMPONENTS

5 Reusable Blades



GlideScope Video Monitor

TFT Color, VGA 640 x 480 px
Monitor: 6.4 in
Height: 183.3 mm
Width: 223 mm
Depth: 79 mm



Video Cable

Length: 2190 ± 55 mm
Diameter: 5.4 mm
Note: Cable shortened for ...

Video Cable



GlideRite Rigid Stylet

Recommended to facilitate successful intubations.
The angle of the stylet complements the unique angle of the LoPro instruments. For use with 6.0 mm or larger endotracheal tubes.



Mobile Stand

Representative Mobile Stand

